



## CT-100

The CT-100 offers all functions needed in a combined table with compact design and easy-to-operate system.



### Technical parameters

1. rocker arm on lap angle:  $\pm 30^\circ$
2. table size (mm): 670 × 400
3. rotary table: 90°
4. movements seating :  
500-700mm
5. seat-load-bearing: 150kg

### Features

- Full stainless steel kick plate
- Overhead lamp with jointed arm
- Chair elevation controls
- Counterbalanced vision tester arm
- Full autoswitch functions

### Electrical parameters

1. input voltage:  
110V/220V AC 50Hz ~ 60Hz
2. electric power: 150W 220V AC
3. lighting: 16W 110V/220V AV
4. Fuse: 5A 250V AC

### Installation size

1. size:  
990mm × 585mm × 2110mm
2. the whole weight: 95Kg