

CT-150

Combined Table

- The new CT-150 Combined Table combines the modern design with the unique ability to effortlessly control the entire exam room environment.

Features:

- Power for ophthalmoscope and retinoscope voltage selectable, 2.5V, ~12V,
- The comfortable and beautiful chair with electric elevator, stable and low noise.
- Easy operation with a well positioned control panel.
- Unique design of arc shaped worktable, creative and useful.

Specifications:

- Arm tilting angle: $\pm 30^\circ$;
- Table rotate angle: 90° ;
- Slide-board: 500mm;
- Chair loading: 150KGS;
- Chair height: 500-700mm;
- Weight: 200Kgs;
- Power: 120W, 220V, 50Hz/110V, 60Hz;
- Illumination bulb: 15W;
- Table size: 840mm \times 405mm



Ray Vision International Corp.

Tel: 86-10-87765481 Fax: 86-10-87765648

Email: sales@ray-vision.com www.ray-vision.com