

# CT-L470 Ophthalmic Unit

- The rotation angle is 360 degrees. It enables multiple devices operating simultaneously.
- The newly added liner switch is an innovation for Y shape table unit, safer and more convenient.
- It offers chargers for ophthalmoscope and retinoscope.

## SPECIFICATIONS

Unit size: 1200\*1055\*1850mm

Table Size: 420\*510mm

Table rotation range: 360°

Swing arm oscillation angle:  $\pm 35^\circ$

Chair up and down range: 540-690mm

Chair loading ability: 200-250kg

Adjustable backrest range: 160°

Electric Voltage: AC220V/110V 50/60HZ

Weight: 200kg

