



Ray Vision™

Ray Vision International Corporation

Add: Rm. 103, Bldg. 3, No. 11,
Dewai Street, Beijing 100088, China

Tel: +86-10-51732722 / 23

Fax: +86-10-51732721

Instrument Table

TB-V300

The TB-V300 Instrument Table features sturdy construction with a new ergonomic design and smooth caster movements with locks.



Height: 640-860mm adjustable

Table top: 40cm×60cm

Bottom: 45cm×52cm

Weight: 25 Kgs

Maximum Load: 50Kgs

Input Voltage: AC220/110, 50/60Hz

Noise lower than 55db in elevation

Rayvision