

SL-R3

Digital Slit Lamp Microscope

■ Features:

Classical Haag-Streit type slit lamp Excellent optics with multi-layered, multi-coated system Three convenient magnification ratios of 10X, 16X and 25X selectable with Galilean-type drum Single-handed elevation system and easy-to-reach tilt mechanism Optional external illumination system.

■ Specifications

■ MICROSCOPE

Type: Galilean-Type
Magnification Changing System: Three position rotating drum
Eyepieces: 12.5×
Angle between eyepieces: 13°
PD Range: 55mm~78mm
Dioptric Adjustment: ±6D

■ SLIT ILLUMINATION

Slit Width: 0~14mm, continuously variable
(at 14mm, slit becomes a circle)
Slit Length: 1~14mm, continuously variable
Lamp: 6V/20W Halogen Lamp
Slit Rotation: 0~180°, continuously adjustable both vertical and horizontal
Slit inclination: 4 step: 5°, 10°, 15°, 20°
Filters: Thermal safety, UV, Red-free, Cobalt Blue
Aperture diameters: 14mm, 10mm, 5mm, 3mm, 2mm, 1mm, 0.2mm



Ray Vision International Corp.

Tel: 86-10-51029818 Fax: 86-10-51029817
Email: sales@ray-vision.com www.ray-vision.com