



# SL-R3L

## Slit Lamp Microscope

### ■ Features:

- Classical Galilean-Type slit lamp.
- Excellent optics with multi-layered, multi-coated system.
- Three convenient magnification ratios of 10X, 16X and 25X selectable with Galilean-type drum.
- Single-handed elevation system and easy-to-reach tilt mechanism.
- Optional external illumination system.

### ■ Specifications

#### Microscope

Type: Galilean-Type

Magnification Changing System: Three position rotating drum

Eyepieces: 12.5x

Angle between eyepieces: 13°

PD Range: 52mm~78mm

Dioptic Adjustment: ±6D

#### Slit Illumination

Slit Width: 0~14mm, continuously variable (at 14mm, slit becomes a circle)

Slit Length: 1~14mm, continuously variable

Lamp: LED bulb

Slit Rotation: 0~180°, continuously adjustable both vertical and horizontal

Slit inclination: 4 step: 5°, 10°, 15°, 20°

Filters: Thermal safety, UV, Red-free, Cobalt Blue

Aperture diameters: 14mm, 10mm, 5mm, 3mm, 2mm, 1mm, 0.2mm

