



SL-R5L

Slit Lamp Microscope

The NEW 3rd Generation of professional slit lamp, with complete function and advanced design for ophthalmologic examination.

■ Specifications

Microscope

Type: Galilean-Type

Magnification Changing System: Five position rotating drum

Eyepieces: 12.5x

Total Magnifications: 6x, 10x, 16x, 25x, 40x

Angle between eyepieces: 13°

PD Range: 52mm~78mm

Dioptic Adjustment: $\pm 6D$

Real Fields of View:

40x (D=5.5mm) , 25x (D=8.5mm), 16x (D=13.5mm),

10x (D=22mm), 6x (D=34.7mm)

Slit Illumination

Slit Width: 0~14mm, continuously variable (at 14mm, slit becomes a circle)

Slit Length: 1~14mm, continuously variable

Lamp: LED bulb

Slit Rotation: 0~180°, continuously adjustable both vertical and horizontal

