

## SL-RVP *Slit Lamp*

### Features:

- Rechargeable and cordless
- Compact, lightweight and ergonomic design
- Easy selection of slit width
- Extra-bright LED lamp light source
- Safe and comfortable for patient
- Fits into lightweight case



### Specifications:

Magnification:	10×, 16× eyepieces included
Diopter Adjustment::	±5D
Diameter of light:	Φ 0.1, Φ 1, Φ 3, Φ 10mm
Pupil distance:	45-70mm
Working distance:	78mm
Slit width:	0mm~ 10mm continuously adjustable
Slit height:	Φ 0.2, Φ 0.3, Φ 5, Φ 10mm
Filter:	Heat absorption, Red-free, Cobalt blue
Illumination rotator angle:	Horizontal ±30°
Illumination:	High luminance white LED