

SL-V2

The SL-V2 combines traditional design with excellent optics to exceed your expectations.

Features

- Traditional Haag-Streit type slit lamp
- Designed for routine practice
- Rapid magnification change from 10X, 16X
- Halogen illumination standard
- Optional Applanation Tonometer for eye care professionals



Specifications:

Microscope Type	Converging stereoscope
Magnifications	Two-step magnifications
Eyepiece	10x, 16x
Magnification (Field of View)	10x (φ18mm), 16x (φ14.5mm), 16x (φ11.25mm), 25.6x (φ9mm)
Range of Pupil Distance(mm)	55mm~82mm (10x), 51mm~78mm(16x)
Adjustment of Diopter	±8D(10x), ±10D(16x)
Slit Projection	2/3x
Slit Width	0mm~9mm adjustable(Round slit at φ9mm)
Slit Height	1mm~8mm adjustable
Aperture Diameter	φ9mm, φ8mm, φ5mm, φ3mm, φ2mm, φ1mm, φ0.2mm
Slit Angle	0°~180°adjustable
Filter	Heat absorption, grey, redfree, cobalt blue
Illumination Bulb	12v30w halogen bulb
Herby Lens	-58.7D
Fixation Bulb	Red LED