

# SL-R5

## Slit Lamp Microscope

- The NEW 3rd Generation of professional slit lamp SL-R5, with complete function and advanced design for ophthalmologic examination.



# SL-R5

## Slit Lamp Microscope

### ■ Features:

- Classical Haag-Streit type slit lamp
- Outstanding Galilean optics with moistureproof, mildewproof and anti-reflection lens
- 0~14mm slit width
- Drum change magnifications 6X, 10X, 16X, 25X and 40X
- Single-handed elevation system and easy-to-reach tilt mechanism
- Optional external illumination system
- Versatility and ease of operation

Specifications	
Type	Galilean-Type
Total Magnifications	6X, 10×, 16×, 25×, 40X (professional configuration)
PD Range	55mm~78mm
Dioptric Adjustment	±6D
Slit Width	0~14mm, continuously variable
Slit Length	1~14mm, continuously variable
Lamp	6V/20W Halogen Lamp
Slit Rotation	0~180° , continuously adjustable both vertical and horizontal
Slit inclination	° 4 step: 5°、10°、15°、20°
Filters	Thermal safety, UV, Red-free, Cobalt Blue
Aperture diameters	∅14mm, ∅10mm, ∅5mm, ∅3mm, ∅2mm, ∅1mm, ∅0.2mm

