

# TL-E80

TL-E80 **Perimeter Trial Lens Set** with leatheroid case contains 2 x 28 pieces Spheres (concave and convex), and 12 pieces of Minus Cylinder. The lenses with silver metal rim is ideal for visual field test.



## SPECIFICATIONS:

- **Spheres**

Lenses with silver metal rim, black leatheroid case

28 pieces each of concave and convex

0.25D to 3.00D in 0.25 steps

3.00D to 7.00D in 0.50 steps

7.00D to 14.00D in 1.00 steps

14.00D to 16.00D in 2.00 steps

- **Minus Cylinder**

12 pieces

0.25D, 0.50D, 0.75D, 1.00D, 1.50D, 2.00D,

2.50D, 3.00D, 3.50D, 4.00D, 4.50D, 5.00D

Ray Vision International Corp.

Tel: 86-10-87765481 Fax: 86-10-87765648

Email: [sales@ray-vision.com](mailto:sales@ray-vision.com) [www.ray-vision.com](http://www.ray-vision.com)