

# TRIAL LENS SET TL-G40

The full diameter trial lens set TL-G40 with color metal rim is very exquisite. The rims and handles of the plus sphere are plated black, and they can be clearly distinguished from the Minus lenses (plated red), even in a dim environment such as an exam room.



## ■ Features:

- ▶ Lenses with color metal rim
- ▶ 266/272 pieces lens as standard with extra 2 pieces of cross cylinder
- ▶ Convenient carrying case made of aluminum
- ▶ Excluding trial frame

## ■ Specifications

### ▶ Spheres

- 40 pairs each of concave and convex
- 0.25D to 6.00D in 0.25 steps
- 6.50D to 10.00D in 0.50 steps
- 11.00D to 15.00D in 1.00 steps
- 16.00D to 20.00D in 2.00 steps

### ▶ Cylinders

- 20 pairs each of concave and convex
- 0.25D to 6.00D in 0.25 steps
- 4.50D to 6.00D in 0.50 steps

### ▶ Prisms

- 12 pieces
- 0.5(2) 1.0 to 10.0 in 1.0 steps

### ▶ Accessories

- 14 pieces
- RF GF BL PL(2) MR(2) FL CL(2) PH(2) SS
- PF Cross cylinder x 0.25 x0.50

