



VISION TESTER VT-350

The VT-350 is a practical vision tester. It can provide quick and accurate refractions very easily. With synchronised cross cylinders, it automatically rotates in conjunction with the set cylinder axis.



■ Specifications

- **Sphere Power Range:** +16.75D ~ -19.00D by 0.25D step
- **Cylinder Power Range:**
0.00D ~ -6.00D
0.00D ~ -8.00D by 0.25D step
(when applied with auxiliary lens-2.00D)
- **Cylinder axis Range:** 0° ~ 180° by 5° step
- **Cross cylinder:** ±0.25D
- **Rotary prism Range:** 0△ ~ 20△ by 1△ step
- **PD adjustment:** 50 ~ 75mm
- **Convergence adjustment:** ∞, 380mm (when PD=64mm)
- **Forehead rest adjustment:** 16mm
- **Vertex distance:** 12mm
- **Standard accessorial lens:** circinal lens, double even lens
- **Standard accessories:** near point card, reading rod, dust cover, face mat, lens accessorial instrument.
- **Size:** 318(L)X96(W)X292(H) (mm)
- **Weight:** ≈4.5kg



Ray Vision International Corp.

Tel: 86-10-87765481 Fax: 86-10-87765648

Email: sales@ray-vision.com www.ray-vision.com