

CP-300

Chart Projector

CP -300 Chart Projector provides improved sharp and clear image projection and fast chart rotation with high-intensity white LED illumination system and instant random access functionality.

Features

- New maintenance-free LED illumination system produces the brighter, more consistent high-resolution image.
- 30 popular optotypes included.
- Masking options allow you to isolate vertical and horizontal lines, or individual characters.
- Fully programmable chart selections.
- Easy setup and maintenance.

Specifications:

- Variable focus projection distance: 1.5m~6m
- Projection magnification: 30 X (in 5m refraction)
- Projection size: 330 X 270mm, ϕ 300mm
- Number of charts: 30
- Chart change-over: 1 frame / 0.03 sec.
- Number of masks: Open 1, Horizontal line 5, Vertical line 8, Single isolation 21, R&G 1.
- Mask change-over: 1 frame / 0.03 sec.
- Program step: 2 type, Max. 30 steps available.
- Light source: LED illumination system.
- Auto-off function: After 10 minutes idle time.
- Power source: AC 220V, 50Hz or 110V, 60Hz
- Power consumption: 80W
- Dimensions: 300(L) X 230(W) X 240(H) mm
- Net weight: 6.0 kg
- Standard Accessories: Remote control, Polarized metal screen, Halogen lamp, Polarized glasses, Fuses(2), Batteries(2)

