

# Elite™

## Digital Visual Acuity Chart




### *Any* SCREEN IS YOUR VISUAL ACUITY CHART !

*High definition, Contrast and sharp display.  
Lifetime of free software upgrades*



- Fits any LCD monitor or TV set.
- High definition, up to 1920 X 1080 (HD 1080P)
- Lifetime of free software upgrades through Wi-Fi.
- Plug and play.
- Hardware and system entirely integrated.
- Multi-functional chart options.
- RF wireless remote, fluorescent keypad.



# Elite™ Digital Visual Acuity Chart



- Smart calibration with easy setup.
- HDMI/VGA supported.
- Fluorescent keypad makes testing in a dark room much easier.
- No direct aiming required, 2.4G RF remote control is omnidirectional.
- Flexible working distance, faster response time, and reliable performance.
- Compact size /w 5G system storage.
- Video, audio, and gallery display.
- Wi-Fi and Bluetooth compatible.
- Continuous testing distance, from 3' to 22'.
- Random sequencing of optotypes.
- Multilingual system along with Imperial and Metric standard.
- Mirror function available for all charts.



Ray Vision International Corp.  
 Tel.: 86-10-51029818  
 Email:sales@ray-vision.com  
 www.ray-vision.com