

# TRIAL LENS SET TL-L15

With the large diameter (26mm) plastic rim, TL-L15 trial lens set adopts the innovative design, fixing tightly the colorful plastic rim with the sliver ring hoop to prevent the plastic rim from breaking, and this trait makes the TL-L15 a typical one with durable usage and ease mobility.



## ■ Features:

- ▶ Lenses with plastic rim
- ▶ 266/272 pieces lens as standard with extra 2 pieces of cross cylinder
- ▶ Convenient carrying case made of wood
- ▶ Excluding trial frame

## ■ Specifications

### ▶ Spheres

- 40 pairs each of concave and convex
- 0.25D to 6.00D in 0.25 steps
- 6.50D to 10.00D in 0.50 steps
- 11.00D to 15.00D in 1.00 steps
- 16.00D to 20.00D in 2.00 steps

### ▶ Cylinders

- 20 pairs each of concave and convex
- 0.25D to 6.00D in 0.25 steps
- 4.50D to 6.00D in 0.50 steps

### ▶ Prisms

- 12 pieces
- 0.5(2) 1.0 to 10.0 in 1.0 steps

### ▶ Accessories

- 14 pieces
- RF GF BL PL(2) MR(2) FL CL(2)PH (2) SS
- PF Cross cylinder x 0.25 x0.50

