

TRIAL LENS SET TL-V20

The large diameter (26mm) plastic rim TL-V20 trial lens set adopts the innovative design. The colorful plastic rim is fixed tightly with the **silver ring hoop** to prevent the plastic rim from breaking, and this trait makes the TL-V20 a typical one with durable usage and ease mobility.



■ Features:

- ▶ Lenses with plastic and metal combi rim
- ▶ 266/272 pieces lens as standard with extra 2 pieces of cross cylinder
- ▶ Convenient carrying case made of aluminum
- ▶ Excluding trial frame

■ Specifications

▶ Spheres

- 40 pairs each of concave and convex
- 0.25D to 6.00D in 0.25 steps
- 6.50D to 10.00D in 0.50 steps
- 11.00D to 15.00D in 1.00 steps
- 16.00D to 20.00D in 2.00 steps

▶ Cylinders

- 20 pairs each of concave and convex
- 0.25D to 6.00D in 0.25 steps
- 4.50D to 6.00D in 0.50 steps

▶ Prisms

- 12 pieces
- 0.5(2) 1.0 to 10.0 in 1.0 steps

▶ Accessories

- 14 pieces
- RF GF BL PL(2) MR(2) FL CL(2) PH(2) SS
- PF Cross cylinder x 0.25 x 0.50

