

SL-V5 *Slit Lamp*

Features:

- Traditional Haag-Streit type slit lamp
- Advanced converging stereoscopic microscope
- Five-step magnification range from 6X to 40X
- Special-layered, multi-coated to transmit light more efficiently
- Available with a variety of optional accessories

Specifications:

Microscope Type:	Galilean stereoscopic
Magnifications:	Five-step drum magnifications
Eyepiece:	12.5x
Magnification Ratio:	6x (φ33) , 10x (φ22.5) , 16x (φ14) , 25x (φ8.8) , 40x (φ5.5)
Range of Pupil Distance (mm) :	55~75
Adjustment of Diopter:	-5D~+3D
Slit Inclination:	5°, 10°, 15°, 20°
Aperture Diameter (mm) :	φ9, φ8, φ5, φ3, φ2, φ1, φ0.2
Filter:	Heat absorption, grey, redfree, cobalt blue
Illumination Bulb:	12V/30W halogen bulb
Optional Accessories:	Applanation Tonometer, Measuring Eyepiece, Trihedron lens, Gonioscope, Funduscope, Aspheric Lens, Teaching Equipment, Hrubby Lens(-58.7D)

